
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

c-Ret (Long Isoform) Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01082

第三版

描述:The ret proto-oncogene products (proto-Ret protein) are expressed as 150kDa and 170kDa glycoproteins in neuroblastoma cells and as 150kDa and 190kDa glycoproteins in leukemia cells.

抗原:Synthetic peptide for C-Terminal of Human c-Ret Long Isoform (ANWMLSPSAAKLMDTFDS)

配方:Lyophilized product from 1% BSA in PBS containing 0.05%NaN₃

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性:Long Isoform specific