

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

C4.4A (301) Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01086

第三版

**描述:**C4.4A is known to be present in the placenta, skin, esophagus and peripheral mononuclear cells in humans. It is a GPI (glycosylphosphatidylinoditol) binding protein which has a similar structure to urokinase receptor (uPAR) and belongs to Ly-6 family as uPAR. While most of Ly-6 family molecules are consisted of single domain, uPAR is of three cysteine-rich domains. C4.4A is a membrane protein which binds to cell membrane in GPI binding site, and is consisted of two cysteine-rich domains and a cysteine lacking third domain. C4.4A was isolated from metastatic rat pancreatic cancer cell lines in 1989, and its human homolog was isolated in 2001. In human malignant tumors, it has been reported that C4.4A is expressed in malignant melanoma, urothelial, lung, esophageal and colon cancer. It is considered that C4.4A mediates invasion and metastasis of cancer cells by binding with laminin and degrading extracellular matrix in the progress of cancer. Therefore, immunological detection of C4.4A in cancerous tissue is expected to be an indicator of tumor metastasis.

**抗原:**Synthetic peptide of a part of Human C4.4A (AGHQDRSNSGQYPAKG).

**配方:**Lyophilized product from PBS containing 1% BSA and 0.05% NaN<sub>3</sub>

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5 μg/ml.

**特异性:**