

---

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

---

CD109 (1383) Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01090

---

第三版

**描述:**CD109 is a glycosyl-phosphatidyl-inositol (GPI) - anchored glycoprotein about 180 - 190 kDa. It is shown that CD109 is expressed in vascular endothelial cells, some epithelial cells, activated T-cells and platelets, and subset of CD34+ megakaryocytic leukemia cells.

**抗原:**Synthetic peptide of a part of Human CD109

**配方:**Lyophilized product from 1% BSA in PBS containing 0.05% NaN<sub>3</sub>

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5 μg/ml.

**特异性:**