
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

CDCrel-1 (C354) Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01092

第三版

描述: Parkinson's Disease (PD) is a relatively common neurodegenerative disorder, which is characterized by the loss of midbrain dopamine (DA) neurons and the presence of Lewy bodies, proteinaceous cytoplasmic inclusions that are abundantly enriched in ubiquitin.

抗原: Synthetic peptide of part of C terminal of human CDCrel-1 (RMQEMLQRMKQQMQDQ)

配方: Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μg/ml.

特异性: Cross-reacts with Mouse and Rat