

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

CIS3/SOCS-3 (C204) Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01096

第三版

描述:The Janus family of protein tyrosine kinases (JAKs) and STAT transcription factors regulate cellular processes involved in cell growth, differentiation, and transformation through their association with cytokine receptors.

抗原:Synthetic peptide for a C-terminal part of Human CIS3/SOCS-3.

配方:Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃.

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性:Antigen is used in same amino acid sequences with human, mouse and rat. Cross reacts with mouse, but not tested cross reactivity with rat.