

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Claudin-12 (C) Anti-Mouse Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01100

第三版

**描述:**The tight junction is identified as a belt-like region in which two lipid-apposing membranes lie close together (tight junction strands). Tight junction strands of the adjacent cells form tightly connected pairs.

**抗原:**Synthetic peptide of part of C terminal of Mouse Claudin-12 (GMHTYSQPYSSRSLSAI )

**配方:**Lyophilized product from 1% BSA in PBS containing 0.05% NaN<sub>3</sub>

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5 μg/ml.

**特异性:**Claudin-12 specific. Not cross reacts with Claudin-1, -2, -3, -4, -5, -6, -7, -8 and -15