
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Claudin-15 (C) Anti-Mouse Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01102

第三版

描述:In human, the claudin superfamily consists of at least 18 members, which are involved on paracellular transport as structural and functional components of tight junction. Claudins are directly associated with ZO-1, 2 and 3 and indirectly with AF-6 and cingulin.

抗原:Synthetic peptide of part of C terminal of Mouse Claudin-15 (TSDESDISFGKYGKNA)

配方:Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性:Claudin-15 specific. Not cross reacts with Claudin-1, -2, -3, -4, -5, -6, -7, -8, and -12