
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Cre recombinase Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01118

第三版

描述: Cre is a 38 kDa recombinase protein from bacteriophage P1 which mediates intramolecular and intermolecular site specific recombination between loxP sites. The role of Cre is to resolve dimers of P1 that arise after replication in order to allow partitioning of the two P1 molecules at cell division.

抗原: Synthetic peptide for middle portion (MLHRRSGLPRPSDSNAV) of Cre recombinase.

配方: Lyophilized product from 1% BSA in PBS containing 0.05%NaN3

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Confirmed by western blotting using Cre Cos/TNE sup.