

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

CRS (CysteinyI tRNA synthetase)

多克隆抗体

产品货号: PA01120

第三版

描述: CRS is a class 1 aminoacyl-tRNA synthetase, cysteinyI tRNA synthetase. Each of the twenty aminoacyl-tRNA synthetase catalyzes the aminoacylation of a specific tRNA or tRNA isoaccepting family with the cognate amino acid.

抗原: Recombinant N-term 6x His-human CRS.

配方: Anti-CRS rabbit serum containing 0.05% sodium azide for preservative.

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: This antibody reacts specifically with 85kDa cysteinyI tRNA synthetase protein.