

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Daple Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01123

第三版

**描述:** Classical (  $\beta$ -catenin dependent) and non-classical (  $\beta$ -catenin non-dependent) Wnt signal pathway are signal pathways that regulate forming and maintaining of embryonic growth or cell structures of organs and amplifying of cells. It is also involved in various diseases including cancers. Especially, it has been well known that the non-classical Wnt signal pathway regulates polarity determination of tissues and migration of cells. It has been revealed that scaffold molecule “ Dishevelled (Dvl) ” and low molecular weight GTPase “ Rho molecular group (Rac, Rho) ” that activates downstream have essential roles in this pathway. Daple is a molecule that was identified by A. Kikuchi et al, Osaka University and it shows that it has a homological sequence with Girdin or Gipie. It is suggested that Daple binds to Dvl and controls activation of Rac depended on Wnt5a stimulus via its mutual interaction. It was discovered that Daple is essential for cell migration and restructuring of actin structure and its molecular mechanism has also an important role for healing of wound of skin by Daple knockout mouse analysis.

**抗原:** Synthetic peptide of Daple (EPGGDPQTVWYEYG).

**配方:** Lyophilized product in PBS containing 1 % BSA and 0.05 % NaN<sub>3</sub>

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$ g/ml.

**特异性:** Cross-react with Human Daple.