

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

dbpA Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01125

第三版

描述:DbpA (DNA Binding protein A) belongs to the Y-box binding protein family that contain a DNA binding domain, called the cold shock domain, of about 80-amino-acid residues. DbpA was identified as the protein binding to the epidermal growth factor receptor enhancer or c-erbB-2 ptomoter. These Y-box binding proteins are reported to have maultiple function, such as the regulations of transription and translation. It was reported that dbpA was a candidate molecule to accelerate the process of the inflammation-induced hepatocarcinogenesis (1), the association with the advanced stages of hepatocellular carcinoma and its nuclear localization as a marker of poor prognosis (2).

抗原:Synthetic peptide of the C terminal of Human dbpA (TENPAPPTQQSSAE)

配方:Lyophilized product from 1% BSA in PBS containing 0.05% NaN3

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: