

---

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

---

DNMT-1 (DNA Methyltransferase-1) Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01127

---

第三版

**描述:**The DNMT-1 is one of the enzymes responsible for methylation. Methylation of DNA at cytosine residues plays an important role in regulation of gene expression, genomic imprinting and is essential for mammalian development.

**抗原:**Synthetic peptide of C terminus of human DNMT-1 (MLAKARESASAKIKE)

**配方:**Lyophilized product from 1% BSA in PBS containing 0.05% NaN<sub>3</sub>

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5 μg/ml.

**特异性:**Confirmed by western blotting using HeLa, SY and HT-29 cells