
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Dock7 Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01132

第三版

描述: Small GTPase of Rho family is one of signaling factors which control alteration of cellular morphology, cell division and differentiation, and plays an important role in regulation of generation and organogenesis at an individual level. Its functions are strictly regulated and its dysfunction causes many diseases such as immunological disease or neurodegenerative disease.

抗原: Synthetic peptide of a part of Human Dock7 (KELFALHPSPDEEE)

配方: Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μg/ml.

特异性: Cross-reacts with rat.