



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Drebrin A (DAS2) Rabbit IgG Affinity Purity

多克隆抗体

产品货号: PA01136

第三版

描述:Drebrin (developmentally regulated brain protein), isolated and identified from chicken brain as a protein associated with brain development, is conserved over species, and there are two types of isoforms, drebrin A and drebrin E in mammalian such as rat and human.

抗原:Synthetic peptide of the part of Drebrin A (FIKASDSGPSSS)

配方:Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性:Reacts with rat and mouse Drebrin A.