



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Endothelin-1 Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01140

第三版

**描述:**Endothelins are 21 amino acid peptides that are produced ubiquitously by vascular endothelial cells, smooth muscle cells, and other cells in different organs. Endothelins are secreted as big-endothelins that are converted to active proteins by the endothelin-converting enzyme.

**抗原:**Synthetic peptide of a part of ET-1 (A common parts of Bovine, Canine, Human, Mouse, Porcine and Rat)

**配方:**Lyophilized product from 1% BSA in PBS containing 0.05%NaN<sub>3</sub>

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5 μg/ml.

**特异性:**