

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse Phospho-STK11 (Thr189) Polyclonal Antibody

多克隆抗体

产品货号: PA01664

第三版

描述:Essential role in G1 cell cycle arrest. Phosphorylates and activates members of the AMPK-related subfamily of protein kinases. Tumor suppressor.

抗原:Synthesized peptide derived from human LKB1 around the phosphorylation site of Thr189

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Ubiquitously expressed. Strongest expression in testis and fetal liver.