
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Rat Phospho-STAT2 (Tyr690) Polyclonal Antibody

多克隆抗体

产品货号: PA01770

第三版

描述:Signal transducer and activator of transcription that mediates signaling by type I IFNs (IFN-alpha and IFN-beta). Following type I IFN binding to cell surface receptors, Jak kinases (TYK2 and JAK1) are activated, leading to tyrosine phosphorylation of STA

抗原:Synthesized peptide derived from human Stat2 around the phosphorylation site of Tyr690

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性: