

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat Phospho-MERTK/TYRO3 (Tyr749/681) Polyclonal Antibody

多克隆抗体

产品货号: PA01831

第三版

描述:MERTK (MER Proto-Oncogene, Tyrosine Kinase) is a Protein Coding gene. Diseases associated with MERTK include Retinitis Pigmentosa 38 and MERTK-Related Retinitis Pigmentosa. Among its related pathways are GPCR Pathway and Nanog in Mammalian ESC Pluripotency. GO annotations related to this gene include transferase activity, transferring phosphorus-containing groups and protein tyrosine kinase activity. An important paralog of this gene is AXL. TYRO3 (TYRO3 Protein Tyrosine Kinase) is a Protein Coding gene. Diseases associated with TYRO3 include External Ear Carcinoma and Lymphocytic Choriomeningitis. Among its related pathways are GPCR Pathway and Nanog in Mammalian ESC Pluripotency. GO annotations related to this gene include transferase activity, transferring phosphorus-containing groups and protein tyrosine kinase activity. An important paralog of this gene is AXL.

抗原:Synthesized peptide derived from human MerTK/Tyro3 around the phosphorylation site of Tyr749/681

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:MERTK is Not expressed in normal B- and T-lymphocytes but is expressed in numerous neoplastic B- and T-cell lines. Highly expressed in testis, ovary, prostate, lung, and kidney, with lower expression in spleen, small intestine, colon, and liver. TYRO3 is Abundant in the brain and lower levels in other tissues.