

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat Phospho-CTBP1 (Ser422) Polyclonal Antibody

多克隆抗体

产品货号: PA01841

第三版

描述: This gene encodes a protein that binds to the C-terminus of adenovirus E1A proteins. This phosphoprotein is a transcriptional repressor and may play a role during cellular proliferation. This protein and the product of a second closely related gene, CTBP2, can dimerize. Both proteins can also interact with a polycomb group protein complex which participates in regulation of gene expression during development. Alternative splicing of transcripts from this gene results in multiple transcript variants.

抗原: Synthesized peptide derived from human CtBP1 around the phosphorylation site of Ser422

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4)，1% BSA

应用: WB 1 ~ 5 μg/ml.

特异性: Expressed in germinal center B-cells.