

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CDYL2 Polyclonal Antibody

多克隆抗体

产品货号: PA02665

第三版

描述: Chromodomain Y-like protein (CDYL) is a related ubiquitous nuclear protein expressed at moderate levels in most tissues. The gene encoding for the CDYL protein localizes to chromosome 6p25.1. CDYL2 contains 1 chromo domain. The exact function of CDYL2 remains unknown.

抗原: Synthesized peptide derived from the N-terminal region of human CDYL2

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Ubiquitously expressed.