

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CYSLTR1 Polyclonal Antibody

多克隆抗体

产品货号: PA03431

第三版

**描述:** This gene encodes a member of the G-protein coupled receptor 1 family. The encoded protein is a receptor for cysteinyl leukotrienes, and is involved in mediating bronchoconstriction via activation of a phosphatidylinositol-calcium second messenger system. Activation of the encoded receptor results in contraction and proliferation of bronchial smooth muscle cells, eosinophil migration, and damage to the mucus layer in the lung. Upregulation of this gene is associated with asthma and dysregulation may also be implicated in cancer. Alternative splicing results in multiple transcript variants. CYSLTR1 (Cysteinyl Leukotriene Receptor 1) is a Protein Coding gene. Diseases associated with CYSLTR1 include Asthma and Adenoid Hypertrophy. Among its related pathways are RET signaling and Peptide ligand-binding receptors. GO annotations related to this gene include G-protein coupled receptor activity and leukotriene receptor activity. An important paralog of this gene is CYSLTR2.

**抗原:** Synthesized peptide derived from the Internal region of human CysLTR1

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) , 1% BSA

**应用:** WB 1 ~ 5  $\mu$ g/ml.

**特异性:** Widely expressed, with highest levels in spleen and peripheral blood leukocytes. Lower expression in several tissues, such as lung (mostly in smooth muscle bundles and alveolar macrophages), placenta, small intestine, pancreas, colon and heart.