

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CDKN2AIP Polyclonal Antibody

多克隆抗体

产品货号: PA07025

第三版

描述:CDKN2AIP (CDKN2A-interacting protein), also known as CARF, is a 580 amino acid protein that activates p53 via p14 ARF (alternate reading frame)-dependent and independent pathways. CDKN2AIP-dependent activation of p53, a protein that up-regulates growth arrest and apoptosis-related genes in response to stress signals, leads to an enhancement of p53 function. Expression levels of CDKN2AIP and p53 show an inverse relationship that is caused by a negative-feedback control via a proteasome-mediated degradation pathway. CDKN2AIP is expressed ubiquitously across tissue samples and, along with p14 ARF, is localized to the perinucleolar region within the nucleus.

抗原:Recombinant protein of human CDKN2AIP

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) , 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Ubiquitously expressed.