

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human DDX53 Polyclonal Antibody

多克隆抗体

产品货号: PA09270

第三版

描述:This intronless gene encodes a protein which contains several domains found in members of the DEAD-box helicase protein family. Other members of this protein family participate in ATP-dependent RNA unwinding.

抗原:Recombinant protein of human DDX53

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed in testis. Wide expression in various cancer tissues and cancer cell lines.