

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CEL Polyclonal Antibody

多克隆抗体

产品货号: PA09663

第三版

**描述:** The protein encoded by this gene is a glycoprotein secreted from the pancreas into the digestive tract and from the lactating mammary gland into human milk. The physiological role of this protein is in cholesterol and lipid-soluble vitamin ester hydrolysis and absorption. This encoded protein promotes large chylomicron production in the intestine. Also its presence in plasma suggests its interactions with cholesterol and oxidized lipoproteins to modulate the progression of atherosclerosis. In pancreatic tumoral cells, this encoded protein is thought to be sequestered within the Golgi compartment and is probably not secreted.

**抗原:** Synthetic peptide of human CEL

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) , 1% BSA

**应用:** WB 1 ~ 5  $\mu$  g/ml.

**特异性:** Mammary gland and pancreas. Expressed by eosinophils.1 Publication