

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human NOL3 Polyclonal Antibody

多克隆抗体

产品货号: PA12935

第三版

描述: This gene encodes an anti-apoptotic protein that has been shown to down-regulate the enzyme activities of caspase 2, caspase 8 and tumor protein p53. Multiple transcript variants encoding different isoforms have been found for this gene. Arc has been shown to localize to the cytoskeleton of neuronal cells and appears to colocalize with F-Actin, although it may associate with an Actin-associated protein rather than directly with F-Actin. It has been shown that cocaine-stimulated neuronal activity results in increased Arc mRNA levels in striatum.

抗原: Recombinant protein of human NOL3

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Highly expressed in heart and skeletal muscle. Detected at low levels in placenta, liver, kidney and pancreas.