

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human APOA4 Polyclonal Antibody

多克隆抗体

产品货号: PA12940

第三版

描述: Apolipoprotein (apo) A-IV gene contains 3 exons separated by two introns. A sequence polymorphism has been identified in the 3'UTR of the third exon. The primary translation product is a 396-residue preprotein which after proteolytic processing is secreted its primary site of synthesis, the intestine, in association with chylomicron particles. Although its precise function is not known, apo A-IV is a potent activator of lecithin-cholesterol acyltransferase in vitro.

抗原: Recombinant protein of human APOA4

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Synthesized primarily in the intestine and secreted in plasma.