

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ATG12 Polyclonal Antibody

多克隆抗体

产品货号: PA12942

第三版

**描述:**Autophagy is a process of bulk protein degradation in which cytoplasmic components, including organelles, are enclosed in double-membrane structures called autophagosomes and delivered to lysosomes or vacuoles for degradation. ATG12 is the human homolog of a yeast protein involved in autophagy.

**抗原:**Recombinant protein of human ATG12

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**Ubiquitous.