

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human AKAP7 Polyclonal Antibody

多克隆抗体

产品货号: PA12997

第三版

描述: This gene encodes a member of the A-kinase anchoring protein (AKAP) family, a group of functionally related proteins that bind to a regulatory subunit (RII) of cAMP-dependent protein kinase A (PKA) and target the enzyme to specific subcellular compartments. AKAPs have a common RII-binding domain, but contain different targeting motifs responsible for directing PKA to distinct intracellular locations. Three alternatively spliced transcript variants encoding different isoforms have been described. [

抗原: Recombinant protein of human AKAP7

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Expressed in brain, heart, lung, pancreas and skeletal muscle.1 Publication