
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human AGTRAP Polyclonal Antibody

多克隆抗体

产品货号: PA12999

第三版

描述: This gene encodes a transmembrane protein localized to the plasma membrane and perinuclear vesicular structures. The gene product interacts with the angiotensin II type I receptor and negatively regulates angiotensin II signaling. Alternative splicing of this gene generates multiple transcript variants encoding different isoforms.

抗原: Recombinant protein of human AGTRAP

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Ubiquitous but more abundant in kidney, heart, pancreas and thyroid.