

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human AGL5 Polyclonal Antibody

多克隆抗体

产品货号: PA13003

第三版

**描述:** Metalloprotease that mediates tubulin deglutamylation. Specifically catalyzes the deglutamylation of the branching point glutamate side chains generated by post-translational glutamylation in proteins such as tubulins. In contrast, it is not able to act as a long-chain deglutamylase that shortens long polyglutamate chains, a process catalyzed by AGTPBP1/ CCP1, AGL1/CCP4 and AGL4/CCP6.

**抗原:** Recombinant protein of human AGL5

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5 μ g/ml.

**特异性:** Expressed in brain.