

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

, 1-383aa, Human, His tag, E.coli

产品货号: TP01003

第三版

别名: Microtubule-associated protein tau isoform 3, MAPT, DDPAC, FLJ31424, FTDP-17, MAPTL, MGC138549, MSTD, MTBT1, MTBT2, PPN

描述: Tau, also known as microtubule-associated protein tau (MAPT), is a protein that stabilizes microtubules. It is abundant in neurons in the central nervous system and is less common elsewhere. When this protein is defective, and no longer stabilizes microtubules properly, it can result in dementias, such as Alzheimer's disease. Recombinant human Tau protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography.

配方: Liquid. 20mM Tris-HCl buffer (pH8.0) containing 10% glycerol, 0.1M NaCl

分子量: 42.1 kDa (403aa) confirmed by MALDI-TOF (molecular weight on SDS-PAGE will appear higher)

序列:

MGSSHHHHHSSGLVPRGSHMAEPRQEFVEMDHAGTYGLGDRKDQGGYTMHQDQEGD TDAGLKAE EAGIGDTPSL
EDEAAGHVTQARMVSKSKDGTGSDDKKAKGADGKTKIATPRGAAPPQKQGQANATRIPAKTPPAPKTPPSSGEPKSGD
RSGYSSPGSPGTPGSRRTPSLPTPPTREP KKVAVVRTPPKSPSSAKSRLQTAPVMPDLKNVSKIGSTENLKHQPGGGKVQ
IINKKLDLSNVQSKCGSKDN IKHVPGGGSVQIVYKPV DLSKVT SKCGSLGNIH HKPGGGQVEVKSEKLD FKDRVQSKIGSL
DNITHVPGGGNKKIETHK LTFRENAKAKTDHGAEIVYKSPV VSGDTS PRHLSNV SSTGSIDMVDSPQLATLADEV SASLAK
QGL

纯度: > 95% by HPLC

浓度: 1.0 mg/ml (determined by BCA assay)

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.