

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

ABHD12B,1-255aa, Human, His tag, E.coli

产品货号: TP01017

第三版

别名: Abhydrolase domain containing 12B, BEM46L3, c14_5314

描述: ABHD12B belongs to the serine esterase family. Recombinant human ABHD12B protein, fused to His-tag at N-terminus, was expressed in E.coli.

配方: Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 0.4M Urea, 10% glycerol

分子量: 31.0 kDa (278aa)

序列:

MGSSHHHHHHSSGLVPRGSHMGSMGLGIWHTVPSCRGEDAKGKDCCWYEAALRDGNPIIVYLHGSAEHRAASHRLKLV
KVLSDGGGFHVLSVDYRFGFDSTGKPTTEGLTTDAICVYEWTKARSGITPVCLWGHSLGTGVATNAAKVLEEKGCPVDAI
VLEAPFTNMWVASINYPLLKIYRNIPGFLRTLMDALRKDKIIFPNDEENVKFLSSPLLILHGEDDRTVPLEYGKKLYEIARNAYR
NKERVKMFVIFPPGFQHNLLCKSPTLLITVRDFLSKQWS

纯度: > 95% by HPLC

浓度: 1 mg/ml (determined by Bradford assay)

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.