



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

ADRM1, 1-407aa, Human, His tag, E.coli

产品货号: TP01076

第三版

别名: Proteasomal ubiquitin receptor ADRM1, ARM1, GP110, Rpn13.

描述: ADRM1, also known as proteasomal ubiquitin receptor ADRM1, functions as a proteasomal ubiquitin receptor. This protein is an integral plasma membrane protein which promotes cell adhesion. ADRM1 is thought to undergo O-linked glycosylation.

配方: Liquid. 20mM Tris-HCl buffer (pH8.0) containing 20% glycerol, 0.1M NaCl, 1mM DTT

分子量: 44.7 kDa (431aa) confirmed by MALDI-TOF

序列:

MGSSHHHHHHSSGLVPRGSHMGSHMTTSGALFPSLVPGSRGASNKYLVEFRAGKMSLKGTTVTPDKRKGLVYIQQTDDSLIHFCWKDRTSGNVEDDLIIFPDDCEFKRVPQCPSGRVYVLKFKAGSKRLFFWMQEPKTDQDEEHCRKVNEYLNPPMPGALGASGSSGHELSALGGEGGLQSL LGNMSHSQLMQLIGPAGLGGLGGLGALTGPGLASLLGSSGPPGSSSSSSRSQSAAVTSSSTTSSTRATPAPSAPAAAASATSPSPAPSSGNGASTAASPTQPIQLSDLQSILATMNV PAGPAGGQQVDLASVLTPEIMAPILANADVQERLLPYLPSGESLPQT ADEIQNTLTSPQFQQALGMFSAALASGQLGPLMCQFGLPAEAVEAANKGDVEAFKAMQNNAKPEQKEGDTKDKKDEEEDMSLD

纯度: > 95% by HPLC

浓度: 1mg/ml (determined by Bradford assay)

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.