

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

ahsg, 19-345aa, Mouse, His tag, Insect cell

产品货号: TP01088

第三版

别名: Alpha-2-HS-glycoprotein isoform 1, Alpha-2-HS-glycoprotein

描述: Ahsg, also known as Alpha-2-HS-glycoprotein isoform 1, is probably involved in differentiation. Recombinant mouse ahsg, fused to His-tag at C-terminus, was expressed in insect cell and purified by using conventional chromatography techniques.

配方: Liquid. In phosphate buffered saline (pH7.4), 10% glycerol.

分子量: 36.1kDa (333aa) 40-57kDa (SDS-PAGE under reducing conditions.)

序列:

APQGTGLGFRELACDDPEAEQVALLAVDYLNHLLQGFKQVLNQIDKVKVWSRRPFGVVYEMEVDTLETTCHALDPTP  
LANCSVRQLTEHAVEGDCDFHILKQDQGFRVMHTQCHSTPDSAEDVRKLCPRCPLLTPFNDTNVVHTVNTALAAFNTQ  
NNGTYFKLVEISRAQNVPLPVSTLVEFVIAATDCTAKEVTDPAKCNLLAEKQHGFCANLMHNLGGEEVSVACKLFQTQ  
PQPANANAVGPVPTANAALPADPPASVVVGPVVVPRGLSDHRTYHDLRHAFSPVASVESASGETLHSPKVGQPGAAGPV  
SPMCPGRIRHFKIHSHHHHH

纯度: > 95% by HPLC

浓度: 0.25mg/ml (determined by absorbance at 280nm)

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.