

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

AHSP, 1-102 aa, Human, Recombinant, E.coli

产品货号: TP01090

第三版

别名: Alpha-hemoglobin stabilizing protein, ERAF, EDRF.

描述: AHSP (Alpha-hemoglobin stabilizing protein), also known as ERAF (Erythroid associated factor), is an erythroid-specific protein that acts as a chaperone to prevent the aggregation of  $\alpha$ -hemoglobin during normal erythroid cell development.

配方: Liquid. In 20 mM Tris-HCl buffer (pH8.0) containing 10% glycerol

分子量: 11.8 kDa (102aa)

序列:

MALLKANKDLISAGLKEFSVLLNQQVFNDPLVSEEDMVTVVEDWMNFYINYYRQQVTGEPQERDKALQELRQELNTLA  
NPFLAKYRDFLKSHELPSHPPSS

纯度: > 95% by HPLC

浓度: 1 mg/ml (determined by Bradford assay)

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.