

---

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

---

AIMP1 (Aminoacyl tRNA synthetase complex-interacting multifunctional protein 1)

产品货号: TP01096

---

第三版

别名:p43, SCYE1, EMAP-2.

描述:AIMP1, previously known as p43, is an auxiliary component of ARS multi complex and a precursor of cytokine endothelial-monocyte activating polypeptide II (EMAPII). AIMP1 is known as a cytokine working in the control of angiogenesis, inflammation and wound healing. AIMP1 also has glucagon-like hormonal activity.

配方:Phosphate buffered saline, pH 7.4

分子量:35kDa

序列:

纯度:> 95% by HPLC

浓度:0.50mg/ml

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.