

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

AIMP3 (Aminoacyl tRNA synthetase complex-interacting multifunctional protein 3)

产品货号: TP01098

第三版

别名:Aminoacyl tRNA synthetase complex-interacting multifunctional protein 3 Protein.

描述:Component of the multisynthase complex which is comprised of a bifunctional glutamylprolyltRNA synthase, the monospecific isoleucyl, leucyl, glutaminy, methionyl, lysyl, arginyl and aspartyl-tRNA synthases, and three auxiliary proteins, EEF1E1/p18, AIMP2/p38 and AIMP1/p43. p18 is a cofactor of multi-ARS complex and regulates p53 through the interaction with ATM/ATR.

配方:1X PBS(pH7.2) 200mM NaCl 20% Glycerol.

分子量:18kDa

序列:

纯度:> 95% by HPLC

浓度:0.50 mg/ml

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.