

本公司提供的电版本说明书仅供参考，实验请以收到的纸质手册为准。

AMMECR1L, 1-310aa, Human, His tag, E.coli

产品货号: TP01146

第三版

别名:AMMECR1-like protein, AMME chromosomal region gene 1-like, MGC4268

描述:AMMECR1L contains 1 AMMECR1 domain. Recombinant human AMMECR1L protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.

配方:Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 10% glycerol, 1mM DTT

分子量:36.9 kDa (333aa) confirmed by MALDI-TOF

序列:

MGSSHHHHHHSSGLVPRGSHMGSMGKRRCVPPLEPKLAAGCCGVKKPKLSGSGTHSHGNQSTTVPGSSSGPLQNHQH  
VDSSSGRENVSDLTLGPGNSPITRMNPNASGALSPLPRPNGTANTTKNLVVTAEMCCYCFDVLYCHLYGFPQPRLPRTNDP  
YPLFVTWKTGRDKRLRGCIGTFSAMNLHSGLREYTLTSALKDSRFPPLTREELPKLFCVSLLTNFEDASDYLDWEVGVHGI  
RIEFINEKGVKRTATYLPVAKQDWDQIQITIDSLLRKGGFKAPITSEFRKTIKLTRYRSEKVTISYAEYIASRQHCFQNGTLH  
APPLYNHYS

纯度:> 95% by HPLC

浓度:0.5 mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.