



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Annexin A5, 1-320aa, Human, E.coli

产品货号: TP01169

第三版

别名:Annexin V, ANX5, ENX2, PP4.

描述:Annexin A5, also known as Annexin V, belongs to the annexin family of calcium-dependent phospholipid binding proteins some of which have been implicated in membrane-related events along exocytotic and endocytotic pathways.

配方:Liquid. In 20 mM Tris-HCl buffer (pH8.0) containing 1mM DTT 10% glycerol

分子量:35.9 kDa (320aa)

序列:

MAQVLRGTVTDFPGFDERADAETLRKAMKGLGTDEESILTLTSSRNAQRQEISAAFKTLFGRDLLDDLKSELTKGFEKLIV  
ALMKPSRLYDAYELKHALKGAGTNEKVLTEIIASRTPEELRAIKQVYEEYGGSSLEDDVVGDTSGYYQRMLVLLQANRDP  
DAGIDEAQVEQDAQALFQAGELKWTGDEEKFITIFGTRSVSHLRKVFVKYMTISGFQIEETIDRETSGNLEQLLLAVVKSIRS  
IPAYLAETLYYAMKGAGTDDHTLIRVMVSRSEIDLFNIRKEFRKNFATSLYSMIKGDTSGDYKKALLLLCGEDD

纯度:> 95% by HPLC

浓度:1 mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.