

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Annexin A6, 1-673aa, Human, His tag, E.coli

产品货号: TP01170

第三版

别名: ANXA6, ANX6, CBP68.

描述: Annexin A6 belongs to a family of calcium-dependent membrane and phospholipid binding proteins. Although their functions are still not clearly defined, several members of the annexin family have been implicated in membrane-related events along exocytotic and endocytotic pathways.

配方: Liquid. 20mM Tris-HCl buffer (pH8.0) containing 10% glycerol 0.1M NaCl

分子量: 78.0 kDa (693aa)

序列:

```
MGSSHHHHHSSGLVPRGSHMAKPAQGAQYRGSIHDFPGFDPNQDAEALYTAMKGFSGDKEAILDIIITSRSNRQRQEV  
QSYKSLYGKDLIADLKYELTGKFERLIVGLMRPPAYCDAKEIKDAISGIGTDEKCLIEILASRTNEQMHQLVAAYKDAYERD  
LEADIIIGDTSGHFQKMLVLLQGTREDDVVSSEDLVQQDVQDLYEAGELKWGTDEAQFIYILGNRSKQHLRLVFDEYLKT  
TGKPIEASIRGELSGDFEKLMLAVVKCIRSTPEYFAERLFKAMKGLGTRDNTLIRIMVSRSELDMLDIREIFRTKYEKSLYS  
MIKNDTSGEYKKTLLKLSGGDDDAAGQFFPEAAQVAYQMWELSAVARVELKGTVRPANDFNPDADAKALRKAMKGLGTDE  
DTIIDIITHRSNVQRQQIRQTFKSHFGRDLMTDLKSEISGDLARLILGLMMPPAHYDAKQLKKAMEGAGTDEKALIEILAT  
RNAEIRAINAYKEDYHKSLEDALSSDTSGHFRRLISLATGHREEGGENLDQAREDAQVAAEILEIADTPSGDKTSLETRFM  
TILCTRSYPHLRRVFQEFIKMTNYDVEHTIKKEMSGDVRDAFVAIVQSVKNKPLFFADKLYKSMKGAGTDEKTLTRIMVSR  
SEIDLLNIRREFIEKYDKSLHQAIEGDTSGDFLKALLALCGGED
```

纯度: > 95% by HPLC

浓度: 0.5 mg/ml (determined by Bradford assay)

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.