

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

ANXA11, 1-505aa, Human, His tag, E.coli

产品货号: TP01180

第三版

别名:Annexin A11, ANX11, CAP50, Annexin-11, Annexin XI.

描述:ANXA11 (Annexin A11) is a member of the annexin family, a group of calcium-dependent phospholipid-binding proteins. Annexins have unique N-terminal domains and conserved C-terminal domains, which contain the calcium-dependent phospholipid-binding sites.

配方:Liquid. In PBS containing 20% glycerol, 2mM DTT, 100mM NaCl

分子量:56.6kDa (525aa) confirmed by MALDI-TOF

序列:

```
MGSSHHHHHHSSGLVPRGSHMSYPGYPPPPGGYPPAAPGGGPWGGAAAYPPPPSMPPIGLDNVATYAGQFNQDYLSGMA
ANMSGTFGGANMPNLYPGAPGAGYPPVPPGGFGQPPSAQQPVPPYGMYPGGNPPSRMPSYPPYPGAPVPGQPMPPPG
QQPPGAYPGQPPVTYPGQPPVPLPGQQQPVPSYPGYPGSGTVTPAVPPTQFGSRGTITDAPGFDPLRDAEVLKAMKGF
TDEQAIIDCLGSRSNKQRQQILLSFKTAYGKDLIKDLKSELNGFEKTILALMKTPVLFDIYEIKEAIKGVGTDEACLEILASRS
NEHIRELNRAYKAEFKKTLEEAIRSDTSGHFQRLISLSQGNRDESTNVDMSLAQRDAQELYAAGENRLGTDESKFNAVLC
SRSRAHLVAVFNEYQRMTGRDIEKSICREMSGDLEEGMLAVVKCLKNTPAFFAERLNKAMRGAGTKDRTLIRIMVSRSET
DLLDIRSEYKRMYGKSLYHDISGDTSGDYRKILLKICGGND
```

纯度:> 95% by HPLC

浓度:0.25mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.