

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

ANXA13, 1-316aa, Human, His tag, E.coli

产品货号: TP01181

第三版

别名:Annexin A13, ANX13, ISA, Annexin-13.

描述:Annexin A13, also known as ANXA13, is a member of the annexin protein family. Members of the annexin family play a role in cytoskeletal interactions, phospholipase inhibition, regulation of cellular growth, and intracellular signal transduction pathways.

配方:Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 30% glycerol, 0.1M NaCl, 1mM EDTA

分子量:37.5 kDa (336aa), confirmed by MALDI-TOF

序列:

MGSSHHHHHHSSGLVPRGSHMGNRHAKASSPQGFVDVDRDAKKNKACKGMGTNEAAIIIEILSGRTSDERQQIKQKYKA
TYGKELEEVLKSELNFEKTALALLDRPSEYAARQLQKAMKGLGTDESVLIEVLCTRTRNKEIIAIKEAYQRLFDRSLESVDK
GDTSGNLKKILVSLQANRNEGDDVDKDLAQDAKDLVDAGEGRWGTDLAFNEVLAKRSYKQLRATFQAYQILIGKDI
EEAIEEETSGDLQKAYLTLVRCADCEYFAERLYKSMKGAGTDEETLIRIVVTRAEVDLQGIKAKFQEKYQKSLSDMVRS
DTSGDFRKLIVALLH

纯度:> 95% by HPLC

浓度:0.5 mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.