

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

ANXA8L2, 1-327aa, Human, His tag, E.coli

产品货号: TP01185

第三版

别名:Annexin A8-like protein 2, ANXA8, bA145E20.2

描述:Annexin A8-like 2, also known as ANXA8L2, is a member of the annexin family of calcium-dependent phospholipid binding proteins. Annexin family members have been implicated as regulators of such diverse processes as ion flux, endocytosis and exocytosis, and cellular adhesion.

配方:Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 0.1M NaCl, 10% glycerol, 1mM DTT

分子量:39.4 kDa (351aa) confirmed by MALDI-TOF

序列:

MGSSHHHHHHSSGLVPRGSHMGSHMAWWKAWIEQEGVTVKSSSHFNPDPAETLYKAMKGIGTNEQAIIDVLTKRSN  
TQRQQIAKSFKAQFGKDLTETLKSELGKFERLIVALMYPYRYEAKELHDAMKGLGTKEGVIIIEILASRTKNQLREIMKAYE  
EDYGSLEEDIQADTSGYLERILVCLLQGSRDDVSSFVDPALALQDAQDLYAAGENIRGTDEMFKFITLCTRSATHLLRVFEE  
YEKIANKSIEDSIKSETHGSLEEAMLTVVKCTQNLHSYFAERLYYAMKGAGTRDGLIRNIVSRSEIDLNLIKCHFCKMYGK  
TLSSMIMEDTSGDYKNALLSLVGSDP

纯度:> 95% by HPLC

浓度:1 mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.