

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

ARF6, 1-175aa, Human, His tag, E.coli

产品货号: TP01209

第三版

别名:ADP-ribosylation factor 6.

描述:ADP-ribosylation factor 6, also known as ARF6, is a member of the ADP ribosylation factor family of GTP-binding proteins. It is localized to the plasma membrane, and regulates vesicular trafficking, remodelling of membrane lipids, and signaling pathways that lead to actin remodeling.

配方:Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 2mM DTT, 20% glycerol, 0.2mM PMSF

分子量:22.2 kDa (195aa), confirmed by MALDI-TOF

序列:

MGSSHHHHHHSSGLVPRGSHMGKVLKIFGNKEMRILMLGLDAAGKTTILYKCLKLQSVTTIPTVGFNVETVITYKNVKF
NVWDVGGQDKIRPLWRHYTGTQGLIFVDCADRDRIDEARQELHRIINDREMRDAIILIFANKQDLDPDAMPHEIQEK
LGLTRIRDRNWYVQPSCATSGDGLYEGLTWLTSNYKS

纯度:> 95% by HPLC

浓度:1 mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.