

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

ARFIP1, 1-373aa, Human, His tag, E.coli

产品货号: TP01211

第三版

别名: Arfaptin-1, HSU52521

描述: ARFIP1 contains 1 AH domain and is a putative target protein of ADP-ribosylation factor. Recombinant human ARFIP1 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.

配方: Liquid. In Phosphate buffered saline (pH7.4) containing 20% glycerol, 1mM DTT

分子量: 44.1 kDa (396aa) confirmed by MALDI-TOF

序列:

MGSSHHHHHHSSGLVPRGSHMGSMASQAESPKNSSAAEIPVTSNGEVDDSDREHSFNRLDKHSLPSGLGLSETQITSHGFDNTK  
EGVIEAGAFQGSPPAPPLPSVMSPSRVAASRLAQQGS DLIVPAGGQRTQTKSGPVILADEIKNPAMEKLELVRKWSLNTYKC  
TRQIISEKLGRGSRVTDLELEAQIDILRDNKKKYENILKLAQTLSTQLFQMVHTQRQLGDAFADLSLKSLELHEEFGYNADT  
QKLLAKNGETLLGAINFFIASVNTLVNKTIEDTLMTVKQYESARIEYDAYRTDLEELNLGPRDANTLPKIEQSQHFLQAHKE  
KYDKMRNDVSVKLFLEENKVKVLHNQLVLFHNAIAAYFAGNQKQLEQLKQFHILKTPGVDAPSWLEEQ

纯度: > 95% by HPLC

浓度: 0.25 mg/ml (determined by Bradford assay)

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.