

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

ARL2BP, 1-163aa, Human, His tag, E.coli

产品货号: TP01223

第三版

别名:ADP-ribosylation factor-like 2 binding protein, BART, BART1.

描述:ARL2BP is an effector of ADP-ribosylation factor-like 2(ARL2) that is essential for nuclear retention of STAT3. This protein binds to ARL2.GTP with high affinity but does not interact with ARL2.GDP, activated ARF, or RHO proteins.

配方:Liquid. In 20 mM Tris-HCl Buffer (pH 7.5) containing 10% Glycerol

分子量:20.9 kDa (183aa), confirmed by MALDI-TOF

序列:

MGSSHHHHHHSSGLVPRGSHMDALEGESFALSFSSASDAEFDVVGYLEDIIMDDEFQLLQRNFMDKYYLEFEDTEENKLI
YTPIFNEYISLVEKYIEEQLLQRIPEFNMAAFTTTLQHHKDEVAGDIFDMLLTFTDFLAFKEMFLDYRAEKEGRGLDLSSGLV
VTSLCKSSSLPASQNNLRH

纯度:> 95% by HPLC

浓度:1 mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.