

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

ARL3, 1-182aa, Human, His tag, E.coli

产品货号: TP01224

第三版

别名:ADP-ribosylation factor-like protein 3, ARFL3.

描述:ARL3, also known as ARFL3, is a member of the ARF-like protein (ARL) subfamily of small GTPases. Unlike ARFs, ARL3 does not activate the cholera toxin ADP-ribosyltransferase. This protein is ubiquitously expressed.

配方:Liquid. 20mM Tris-HCl buffer (pH8.0) containing 30% glycerol 0.1M NaCl,1mM DTT, 0.1mM PMSF

分子量:23.0 kDa (206aa), confirmed by MALDI-TOF

序列:

MGSSHHHHHHSSGLVPRGSHMGSHMGLLSILRKLKSAPDQEVRIILLGLDNAGKTTLLKQLASEDISHITPTQGFNIKSVQ
SQGFKLNVDIGGQRKIRPYWKNYFENTDILYVIDSADRKRFEETGQELAELEEEKLSCVPVLIFANKQDLLTAAPASEIA
EGLNLHTIRDRVWQIQSCSALTGEGVQDGMNWVCKNVNAKKK

纯度:> 95% by HPLC

浓度:1.0 mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.