

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

ARL4A, 1-200aa, Human, His tag, E.coli

产品货号: TP01225

第三版

别名:ADP-ribosylation factor-like protein 4A, ARL4

描述:ARL4A, also known as ADP-ribosylation factor-like protein 4A, is a member of the ARF-like protein (ARL) subfamily of small GTPases. It is most closely related to ARL6 and ARL7. Unlike ARFs, ARL4 does not activate the cholera toxin ADP-ribosyltransferase.

配方:Liquid. 20mM Tris-HCl buffer (pH8.0) containing 20% glycerol 0.1M NaCl,1mM DTT

分子量:24.7 kDa (220aa), confirmed by MALDI-TOF

序列:

MGSSHHHHHHSSGLVPRGSHMGNGLS DQTSILSNLPSFQSFHIVILGLDCAGKTTVLYRLQFNEFVNTVPTKGFNTEKIKV
TLGNSKTVTFHFWDVGGQEKL RPLWKS YTRCTDGIVFVVDSDVERMEEAKTELHKITRISENQGV PVLIVANKQDLRN
SLSLSEIEKLLAMGELSSSTPWHLQPTCAIIGDGLKEGLEKLHDMIIKRRKMLRQQKKKR

纯度:> 95% by HPLC

浓度:0.5 mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.