

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

ASF1A, 1-204aa, Human, His-tagged, Recombinant, E.coli

产品货号: TP01244

第三版

别名:ASF1 anti-silencing function 1 homolog A, Histone chaperone ASF1A, CCG1-interacting factor A (CIA), CGI-98, HSPC146.

描述:ASF1 anti-silencing function 1 homolog A, also known as ASF1A, belongs to H3/H4 family of histone chaperone proteins. ASF1A protein interacts with histones H3 and H4, and functions together with a chromatin assembly factor during DNA replication and repair.

配方:Liquid. In 20 mM Tris-HCl buffer (pH 8.0) containing 1mM DTT , 10% glycerol

分子量:27 kDa (240aa)

序列:

MRGSHHHHHHGMASMTGGQQMGRDLYDDDDKDRWGSMKVQVNNVVVLDNPSPFYNPFQFEITFECIEDLSEDL
WKIIYVGSASEEYDQVLDSVLVGPVPAGRHMVFVQADAPNPGLIPDADAVGVTVVLITCTYRGQEFIRVGYVNNNEYTE
TELRENPPVKPDFSKLQRNILASNPRVTRFHINWEDNTEKLEDAESSNPNLQSLSTDALPSASKGWSTSENSLNMLESHM
DCM

纯度:> 95% by HPLC

浓度:1 mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.