

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

ASPHD1,96-390aa, Human, His tag, E.coli

产品货号: TP01253

第三版

别名:Aspartate beta-hydroxylase domain containing 1, Aspartate beta-hydroxylase domain-containing 1

描述:ASPHD1 belongs to the aspartyl/asparaginyl beta-hydroxylase family. Recombinant human ASPHD1 protein, fused to His-tag at N-terminus, was expressed in E.coli.

配方:Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 0.4M Urea, 10% glycerol

分子量:33.6 kDa (318aa)

序列:

MGSSHHHHHHSSGLVPRGSHMGSYRLGSQDMQALGAGSRAGGVRGGPVGCSEAGGPPSPGGPGDPGEGPRTEGLVSRRL  
RAYARRYSWAGMGRVRRAAQGGPGPRGPGVGLGIQRPGLLFLPDLPAPFVPRDAQRHDELLESSFPAILRDFGAVSWD  
FSGTTPPPRGWSPPLAPGCYQLLLYQAGRCQPSNCRRCPGAYRALRGLRSFMSANTFGNAGFSVLLPGARLEGRGCGPTNA  
RVRCHLGLKIPPGCELVVGGEPCWAEGHCLLVDDSLHTVAHNGSPEDGPRVVFIVDLWHPNVAGAERQALDFVFAP  
DP

纯度:> 95% by HPLC

浓度:0.5 mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.